

## Erbium 2000 Er:YAG Laser

--For Medicine Skin Resurfacing,  
Smoothing Scars (Dermabrasion).



- Modular Design
- OEM Design and Supply
- Good Beam Quality
- Compact Construction
- High Reliability

### SPECIFICATION:

LASER HEAD	
Laser:	Er:YAG
Wavelength:	2.94 $\mu$ m
Max. Average Power:	15 W
Max. Pulse Energy:	1.5 J
Pulse Repetition Rate:	1 Hz --10 Hz / Adjustable
Spot Size:	6.5 mm
Pulse Duration:	350 $\mu$ s
Dimension (L x W x H):	500 x 140 x 100 mm
POWER SUPPLY:	
Electrical:	220-240VAC--16A--50/60 Hz
Dimension (L x W x H):	400 x 400 x 200 mm



Omega Optics  
10435 Burnet Rd, Suite 108  
Austin, TX 78758  
Tel: 512-996-8833 ext 605, 602  
Fax: 512-996-8833  
Email: [sales@eomgaoptics.com](mailto:sales@eomgaoptics.com)  
[www.eomegaoptics.com/lasers](http://www.eomegaoptics.com/lasers)